

superDimension[®] Marker Delivery Kit (SDMDK1000)

Introducing the newest, most convenient approach to implanting fiducial markers to identify diseased lung tissue for future treatments.



Kit Contents

Quantity (1) Catheter

- 5 French (0.070 inch/1.8 mm) outside diameter
- 0.040 inch inside diameter
- 125 cm length

Quantity (1) Guidewire

- 0.038 inch outside diameter guidewire
- 150 cm length

Do not use with fiducial markers that exceed 0.9 mm outer diameter.

Product Features

Innovative

- First marker delivery tool indicated for use with the superDimension[®] System
- Enables marker placement in soft lung tissue

Convenient

- Both catheter and guidewire included in the same package

Easy to Use

- Pre-marked for use with superDimension Extended Working Channel (EWC)



Ordering Information for superDimension® Marker Delivery Kits

| Part Number | Description | Unit |
|-------------|------------------------------------|-----------|
| SDMDK1000 | superDimension Marker Delivery Kit | 5 per box |

Marker Delivery Kits complement superDimension's full line of superLock™ Fiducial Markers



| Part Number | Description | Unit |
|-------------|---|-----------|
| SDM2 | superLock Two-Band Marker (0.9 mm x 13 mm) | 5 per box |
| SDM3 | superLock Three-Band Marker (0.9 mm x 23 mm) | 5 per box |
| SDC5 | superLock 5 mm Coil Marker (0.9 mm x 5 mm) | 5 per box |
| SDC10 | superLock 10 mm Coil (0.9 mm x 10 mm) | 5 per box |
| SDC20 | superLock 20 mm Coil (0.9 mm x 20 mm) | 5 per box |

To order superDimension Marker Delivery Kits, or learn more about how our products and technology can facilitate the diagnosis and treatment of lung disease, please contact your local Sales Representative at 1.888.586.4767.

superDimension®

161 Cheshire Lane
Suite 100
Minneapolis, MN 55441
763.210.4000
888.586.4767

www.superdimension.com